## Brain Signal Processing: Basics of EEG Analysis

 $Palani\ Thanaraj\ K$  2020-03-30

# Contents

Prerequisites									5		
1	1 Introduction 1.1 How EEG is used in cl	linical setting							 		. 7
2	2 What is Schizophrenia c 2.1 Dataset	condition?							 		9
3		EEG data import & Analysis							11		
	3.1 MNE Library for EEG	Hanalysis							 		. 11
	3.2 Data Import								 		. 11
	3.3 Time series plot of the	EEG signal							 		. 12
	3.4 Bandpass filtering of the	he EEG signal							 		. 13
	3.5 The power spectral den	The power spectral density (PSD)function of the EEG signal							. 15		
	3.6 More Links								 		. 17

4 CONTENTS

# Prerequisites

This tutorial provides basics of EEG analysis with minimal coding using Python. So not much of technical knowledge is required for going through this tutorial. As and when some jargon is encountered it is provided with required information. Feel free to go through this tutorial.

6 CONTENTS

### Chapter 1

### Introduction

EEG is the recording of electrical signals obtained from the scalp region of the human head. It is because of the activation of the pyramidal neurons in the cerebral cortex region of the brain (i.e outer layer of the brain).

#### 1.1 How EEG is used in clinical setting

EEG represents the functional operation of the brain. The EEG signal is also referred to as Brain Waves that are named based on the frequency component in the signal. They are Delta, Theta, Alpha, Beta, and Gamma brain waves. Their frequency component falls in the range of 1-3,3-7,7-13,13-30 and >30 Hz respectively. EEG is used in the management of brain disorders such as Epilepsy, Schizophrenia (SZ) and Alzhemier's disease. In this tutorial, we will discuss about some aspects of Schizophrenia disease and its assessment using EEG signal.

### Chapter 2

## What is Schizophrenia condition?

Schizophrenia is a detrimental condition of the brain which leads to reduction in the cognitive awareness. The patients affected by this disease try to have impairment in assessing day-to-day activities. EEG is an important modality that can detect the SZ occurrence in an individual. The EEG signal shows low randomness compared to a healthy person. Here we will take few EEG signal of normal and SZ patients and analyse them.

#### 2.1 Dataset

The SZ dataset which we are going to analyse here is from "EEG in schizophrenia (RepOD)" from the seminal works of (Olejarczyk and Jernajczyk, 2017a). The dataset comprised 14 patients with paranoid schizophrenia and 14 healthy controls. Data were acquired with the sampling frequency of 250 Hz using the standard 10-20 EEG montage with 19 EEG channels: Fp1, Fp2, F7, F3, Fz, F4, F8, T3, C3, Cz, C4, T4, T5, P3, Pz, P4, T6, O1, O2. The reference electrode was placed between electrodes Fz and Cz.

### Chapter 3

## EEG data import & Analysis

The EEG data is available in the '\*.edf' file extension. We have to download the EEG files and place in the working directory of the Python program such as Jupyter notebook or Spyder. Then we will be using the MNE Python library for the processing of the EEG signal. The following section provides the code block for import the necessary library files for EEG analysis.

#### 3.1 MNE Library for EEG Analysis

```
import mne # MNE library for EEG analysis
import matplotlib.pyplot as plt # Plotting library

plt.style.use('seaborn-bright') # Plot Style
plt.rcParams['figure.dpi'] = 300 # DPI setting
plt.rcParams['figure.figsize'] = (8, 6)
```

### 3.2 Data Import

## EDF file detected

## Setting channel info structure...

This code chunk explains how to import the EEG file into the Python workspace by the command mne.io.read\_raw\_edf. Moreover, it also shows the file info of the EEG file imported into the workspace.

```
# Data import
raw_h = mne.io.read_raw_edf("h01.edf", preload=True)# Import data file into workspace
# h01 specifies the EEG record of healthy control-1

## Extracting EDF parameters from C:\Users\RajuPC\Documents\MyR\MyGitbook\EEG\EEG_analysis\h01.edf...
## EDF file detected
## Setting channel info structure...
## Creating raw.info structure...
## Reading 0 ... 231249 = 0.000 ... 924.996 secs...
raw_s = mne.io.read_raw_edf("s01.edf", preload=True)# Import data file into workspace
# s01 specifies the EEG record of SZ patient-1

# Give the sample rate

## Extracting EDF parameters from C:\Users\RajuPC\Documents\MyR\MyGitbook\EEG\EEG_analysis\s01.edf...
```

## Creating raw.info structure...

```
## Reading 0 ... 211249 = 0.000 ... 844.996 secs...
print('sample rate:', raw_h.info['sfreq'], 'Hz')
# Give the size of the data matrix

## sample rate: 250.0 Hz
print('Data matrix size : %s channels x %s samples' % (raw_h.info['nchan'], len(raw_h.times)))


## Data matrix size : 19 channels x 231250 samples
print('Channel names:', raw_h.info['ch_names'])
```

## Channel names: ['Fp2', 'F8', 'T4', 'T6', '02', 'Fp1', 'F7', 'T3', 'T5', '01', 'F4', 'C4', 'P4', 'F3'

#### 3.3 Time series plot of the EEG signal


This code plots the raw EEG signal using the available function plot. Healthy Individual

```
# Plotting the raw EEG signal
raw_h.plot(start=300,duration=2, block=True, show=False, color='green')
```



```
SZ Patient
```

```
# Plotting the raw EEG signal
raw_s.plot(start=300,duration=2, block=True, show=False, color='red')
```




#### 3.4 Bandpass filtering of the EEG signal

This code applies a FIR filter based on windowing technique to design a bandpass filter with lower cut-off frequency at 0.5Hz and upper cut-off frequency at 100Hz. Healthy Control

```
xfil_h=raw_h.filter(0.5, 100., fir_design='firwin')#Bandpass filtering
```

```
## Filtering raw data in 1 contiguous segment
## Setting up band-pass filter from 0.5 - 1e+02 Hz
##
## FIR filter parameters
## -------
## Designing a one-pass, zero-phase, non-causal bandpass filter:
## - Windowed time-domain design (firwin) method
## - Hamming window with 0.0194 passband ripple and 53 dB stopband attenuation
## - Lower passband edge: 0.50
## - Lower transition bandwidth: 0.50 Hz (-6 dB cutoff frequency: 0.25 Hz)
## - Upper passband edge: 100.00 Hz
## - Upper transition bandwidth: 25.00 Hz (-6 dB cutoff frequency: 112.50 Hz)
## - Filter length: 1651 samples (6.604 sec)


xfil_h.plot(start=300, duration=2, bgcolor='w',show=False, events=None,proj=False,color='green')#plotti
```



#### SZ patient


```
xfil_s=raw_s.filter(0.5, 100., fir_design='firwin')#Bandpass filtering
```

```
## Filtering raw data in 1 contiguous segment
## Setting up band-pass filter from 0.5 - 1e+02 Hz
##
## FIR filter parameters
## -------
## Designing a one-pass, zero-phase, non-causal bandpass filter:
## - Windowed time-domain design (firwin) method
## - Hamming window with 0.0194 passband ripple and 53 dB stopband attenuation
## - Lower passband edge: 0.50
## - Lower transition bandwidth: 0.50 Hz (-6 dB cutoff frequency: 0.25 Hz)
## - Upper passband edge: 100.00 Hz
## - Upper transition bandwidth: 25.00 Hz (-6 dB cutoff frequency: 112.50 Hz)
## - Filter length: 1651 samples (6.604 sec)
xfil_s.plot(start=300, duration=2, bgcolor='w',show=False, events=None,proj=False,color='red')#plotting
```




### 3.5 The power spectral density (PSD)function of the EEG signal

This code performs the Power Spectral Density (PSD) function of the Filtered EEG signal. Healthy Control xfil\_h.plot\_psd(area\_mode='range', tmax=10.0, show=False, average=True, color='green')



xfil\_s.plot\_psd(area\_mode='range', tmax=10.0, show=False, average=True, color='red')

3.6. MORE LINKS



#### 3.6 More Links

More information related to applications of EEG in brain disorders can be found in the following articles:

- (T. Krishnan et al., 2016)
- (Thanaraj and Chitra, 2014)
- (Krishnan and Balasubramanian, 2017)
- (Krishnan et al., 2017)
- (Thanaraj and Parvathavarthini, 2017)
- (Abhinaya and Charanya, 2016)
- (Johannesen et al., 2016)
- (Götz et al., 2017)
- (Thilakvathi et al., 2017)
- (Olejarczyk and Jernajczyk, 2017b)

## Bibliography

- Abhinaya, B. and Charanya, D. (2016). Feature Extraction and Selection of a Combination of Entropy Features for Real-time Epilepsy Detection. *International Journal Of Engineering And Computer Science*, 5(5):16073–16078.
- Götz, T., Stadler, L., Fraunhofer, G., Tomé, A. M., Hausner, H., and Lang, E. W. (2017). A combined cICA-EEMD analysis of EEG recordings from depressed or schizophrenic patients during olfactory stimulation. *Journal of Neural Engineering*, 14(1):016011.
- Johannesen, J. K., Bi, J., Jiang, R., Kenney, J. G., and Chen, C.-M. A. (2016). Machine learning identification of EEG features predicting working memory performance in schizophrenia and healthy adults. *Neuropsychiatric Electrophysiology*, 2(1).
- Krishnan, P. T. and Balasubramanian, P. (2017). Automated EEG seizure detection based on S-Transform. In 2016 IEEE International Conference on Computational Intelligence and Computing Research, ICCIC 2016, pages 1–5. IEEE.
- Krishnan, P. T., Balasubramanian, P., and Krishnan, C. (2017). Evaluation of joint analysis of multiple interictal events for spike extraction and source imaging. *International Journal of Biomedical Engineering and Technology*, 25(1):77.
- Olejarczyk, E. and Jernajczyk, W. (2017a). Eeg in schizophrenia.
- Olejarczyk, E. and Jernajczyk, W. (2017b). Graph-based analysis of brain connectivity in schizophrenia. *PLOS ONE*, 12(11):e0188629.
- T. Krishnan, P., Balasubramanian, P., and Krishnan, C. (2016). Segmentation of Brain Regions by Integrating Meta Heuristic Multilevel Threshold with Markov Random Field. *Current Medical Imaging Reviews*, 12(1):4–12.
- Thanaraj, K. P. and Chitra, K. (2014). Multichannel Feature Extraction and Classification of Epileptic States Using Higher Order Statistics and Complexity Measures. *International Journal of Engineering* {&} Technology (0975-4024), 6(1):102-109.
- Thanaraj, P. and Parvathavarthini, B. (2017). Multichannel interictal spike activity detection using time–frequency entropy measure. Australasian Physical & Engineering Sciences in Medicine, 40(2):413–425.
- Thilakvathi, B., Shenbaga Devi, S., Bhanu, K., and Malaippan, M. (2017). EEG signal complexity analysis for schizophrenia during rest and mental activity. *Biomedical Research*, 28(1):1–9.